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BRAF mutation testing: the different tests

Whilst not an exhaustive list, the table below provides details of the available tests for identifying *BRAF* mutations. Regarding next generation sequencing, although this isn't available as a testing kit, for laboratories with sequencing technologies available this is an alternative to the below.

Test ^a	Method, instrumentation	Regulatory status	Analytical sensitivity LOD	<i>BRAF</i> mutations detected
THxID™- <i>BRAF</i> ¹ (bioMérieux, Inc.)	Real-time PCR, ABI7500 FastDx (open)	FDA ² CE-IVD	5% mutation level for V600E and K	V600E and V600K
cobas® 4800 <i>BRAF</i> V600 ³ (Roche Diagnostics)	Real-time PCR, Cobas 4800 (closed)	FDA ² CE-IVD	5% mutation level for V600E	V600E (V600K and V600D possibly by cross reactivity)
therascreen <i>BRAF</i> RGQ PCR Kit ⁴ (QIAGEN)	Real-time PCR, Rotor-Gene Q HDx 5plex HRM instrument (open)	CE-IVD	1.87% for V600E 4.3% for V600K	V600E, V600E complex, V600D, V600K, V600R
therascreen <i>BRAF</i> Pyro Kit ⁵ (QIAGEN)	Pyrosequencing, PyroMark Q24 System (open)	CE-IVD	Not established	Codon 600: V600A, V600E, V600G, V600M, Codons 464-469: G464E, G464V, G466E, G466V, G469A, G469E, G469V
INFINITI <i>BRAF</i> Assay (AutoGenomics, Inc.) ⁶	INFINITI Analyzer, INFINITI Plus Analyzer	CE-IVD	Not available	<i>BRAF</i> V600A, V600D, V600E, V600KRM
FoundationOne CDx (Foundation Medicine, Inc.) ⁷	NGS, multiple instruments	FDA ²	≤ 2.0% allele fraction for <i>BRAF</i> V600E and K	<i>BRAF</i> V600E, V600K, plus other genes

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